



Implant Cervical Collars: Preserving Crestal Bone

Neil Park, DMD

January 30, 2017



Implant Cervical Collars: Preserving Crestal Bone

ABSTRACT

Roughened implant surfaces have been instrumental in bolstering long-term implant survival rates, but can also promote microbial attachment if exposed to the gingival crevice or oral cavity. The likelihood that some bone loss will occur from surgical trauma or other contributing factors over time is considerable. Consequently, maintaining a machined surface on the implant collar can help to prevent or minimize the percentage of roughened surface exposure after implant placement.

INTRODUCTION

Implant restoration of missing teeth has become increasingly predictable since 2000,¹ when long-term implant survival rates began to exceed 95%.¹⁻⁴ Although a host of internal and external variables can still significantly affect the osseointegration process (Table 1⁵⁻¹⁵), this significant enhancement in implant predictability over the preceding two decades¹ can be attributed to a number of technological improvements in component designs and surgical techniques. Among these, dental implant surfaces are a crucial aspect of the osseointegration process because they directly influence bone and soft tissue attachment to the implant.¹⁶ Numerous studies have reported a positive correlation between implant surface roughness and the overall percentage of bone-to-implant attachment that develops during osseointegration.¹⁷⁻¹⁹ For this reason, most contemporary manufacturers further roughen their dental implant surfaces after machining. However, there is currently no clinical consensus as to what degree of surface roughness is clinically optimal, or even how surface roughness should be measured. For example, while one surface may appear rougher than another surface when average roughness is measured in a two-dimensional (2D) roughness profile (Ra value),¹⁶ the reverse can be true when measuring the same two surfaces using a three-dimensional (3D) roughness profile (Sa value).²⁰ Although some researchers²¹ have strongly asserted that 3D Sa values are more accurate and preferable to measuring implant surface roughness than 2D Ra values, there is still no clinical consensus and both Ra and Sa roughness measurements are used in the industry today.

Despite the recent advancements in implant predictability, crestal bone loss continues to be a matter of concern.²²⁻²⁵ This can be especially problematic following immediate implant placement into fresh extraction sockets, because bone will typically resorb 0.5-1.0 mm in both the vertical and horizontal dimensions during the first year after implant placement due to the post-extraction healing process.^{26, 27} Such bone loss not only harbors a risk of exposing the roughened implant surface to the gingival crevice, but may also directly expose the roughened surface to the oral cavity if there is corresponding gingival recession during the early stage of healing.²⁶ Any degree of roughened surface exposure may also enhance attachment by the 100-200 varieties of intraoral microbiota that inhabit neighboring dentition at any given time.²⁸⁻³⁰ While clinically harmless below a certain count,^{29, 31} unchecked bacterial reproduction results in biofilm formation that can trigger gingival and periodontal diseases in the gingival sulcus regions of both implants and teeth.^{30, 31, 32} While infection is usually manifested as localized gingival inflammation, continued release of inflammatory cytokines, chemokines, and other intracellular mediators from the biofilms can progressively lead to chronic gingivitis, periodontal attachment loss, bone resorption, and eventual tooth loss or implant failure.³⁰⁻³³ This

paper will review issues related to roughened implant surfaces and crestal bone loss as the basis for advocating maintenance of a 1.0 mm machined collar below the prosthetic platform of dental implants.

PERIIMPLANT BONE LOSS

For implants placed into healed extraction sites, approximately 0.5-1.0 mm of periimplant marginal bone loss — often described as *saucerization* — has been observed around implants at their second-stage surgical uncovering. Because surgical insult to the bone tissue during routine implantation procedures is widely assumed to be the cause,^{13, 34-42} such preliminary bone loss is generally accepted as an unavoidable consequence of implant surgery. Disruption of the vascular network during elevation of the mucoperiosteum has also been cited as an additional theoretical cause of this bone loss. The latter theory is not universally accepted, however, because similar saucerization does not appear around natural teeth after soft tissue elevation for osseous surgery,⁴³ and some bone loss still occurs on implants placed via a transmucosal or flapless technique.⁴⁴ Consequently, many implant studies have either failed to report any marginal bone changes,⁴⁵ or only calculated marginal bone loss after definitive implant loading (i.e. at least 3-6 months after implant placement).^{22, 37} As technology improved, early researchers gradually applied new definitions of acceptable periimplant bone loss over time (Table 2^{42, 46-49}).

In addition to the effect of surgical trauma^{13, 34-42} on short-term periimplant bone loss, long-term bone loss and eventual implant failure have also been linked to many other variables, such as bacteria-related infections,^{34, 49-53} chronic inflammation,⁵⁴⁻⁵⁷ and occlusal overloading.^{34, 58-61} In a prospective, multicenter study conducted by the U.S. Department of Veterans Affairs, researchers measured the residual facial plate thickness of approximately 3,000 implant osteotomies and followed patients for three years.⁶² The study found that, as residual facial plates approached 2.0 mm in thickness, there was a corresponding decrease in facial bone loss, bone gain in some patients, and an increase in implant survival rates.⁶² Conversely, as facial plate thickness decreased below the 2.0 mm threshold, there was a corresponding increase in facial bone loss and implant failure.⁶² Implant placement in low-density bone, tobacco use and excessive alcohol consumption have also been associated with increased peri-implant bone loss in the dental literature.^{63, 64}

In the same prospective multicenter study conducted by the U.S. Department of Veterans Affairs, other researchers⁶⁵ evaluated clinician subjective estimates of bone density during osteotomy preparations for 2,839 implants. Results indicated that the greatest bone density was located in the anterior mandible, with diminishing bone density in the posterior mandible, anterior maxilla, and posterior maxilla, respectively.⁶⁵ In an unrelated study, researchers used finite element analysis (FEA) to evaluate bone strains relative to bone properties and simulated periimplant loading conditions.⁶⁶ In cancellous bone models with a lower range of Young's modulus values, FEA results showed that 50% of patients experienced hyperphysiologic periimplant strains in the crestal bone region where short-term periimplant crestal bone loss (saucerization) is commonly reported.⁶⁶ In denser bone with higher ranges of Young's modulus values, only 25% of the cancellous bone models exhibited excessive periimplant bone strains.⁶⁶ These findings suggest that implants placed in low-density bone with a lower range of Young's modulus values characteristic of the anterior and posterior maxillary jaw⁶⁵ may be at greater risk for bone loss than implants placed in other jaw locations with higher Young's modulus values.⁶⁶ It is important to note, however, that low-density bone can also be found in other jaw locations, such as the posterior mandible. Maintaining 1.0 mm of machined surface on the implant collar can help to prevent exposure of the roughened surface in the likely event of some crestal bone loss over time.

IMPLANT CERVICAL MICROARCHITECTURE

Traditional implant dentistry strictly advocated a two-stage, delayed loading technique to prevent implant micromovements⁶⁷⁻⁶⁸ that might result in fibrous tissue encapsulation and clinical failure of the implant. More recently, cellular,⁶⁹⁻⁷¹ animal⁷²⁻⁷⁴ and human⁷⁵⁻⁷⁷ studies have suggested that implant cervical regions with microgrooves,

microthreads, microtextured surfaces or a combination of these features might help to foster bone attachment to the cervical region of the implant surface and bone formation inside the microgrooves or between the microthreads to impede epithelial tissue down-growth. The shapes and dimensions of the microgrooves themselves have also been reported to influence the behavior of epithelial and connective-tissue cells differently *in vivo*.⁷⁸⁻⁸³ Although adequate, prospective, comparative research data on microgrooves and microthreads are still generally lacking, short-term results have been variable in both animal and human models,^{80, 81} while longer term results have suggested that any possible benefit of cervical microthreads may disappear after 5 years in function.⁸² It is currently unknown whether the perceived benefits of cervical microgrooves or microthreads will outweigh the potential hygiene risks of thread exposure should bone recession occur.

Extending roughened implant surfaces to the tops of implant necks has been reported to reduce the amount of periimplant crestal bone loss without adversely affecting soft tissue health.^{79, 83} Implants with roughened, acid-etched (test group) and machined (control group) surfaces were placed in a canine model and allowed to accumulate plaque for 6 months.⁸⁴ There were no significant differences in plaque formation or establishment of inflammatory cell lesions in the periimplant mucosa between the test and control groups.⁸² It should be noted, however, that the acid-etched implants are much smoother than many of the roughened surfaces on the market today.²⁰ Nonetheless, concern about the possibility of increased bacterial attachment to rougher implant surfaces in humans has been expressed in the dental literature.^{79, 83}

RISKS OF ROUGHENED SURFACE EXPOSURE

Concern about possible bacterial colonization of the implant-abutment microgap within the biologic width may be significantly heightened if the implant neck has a roughened surface in close proximity. Although some clinicians have advocated moving the microgap away from the outer circumference of the implant with the use of an abutment that is smaller in diameter than the implant itself (platform switching), the effect of a roughened surface at the microgap remains uncertain. One study⁸⁵ used FEA to evaluate whether platform switching provided any biomechanical advantages. Results showed that, although platform switching helped to shift the stress concentration area away from the cervical bone-implant interface, it increased stress concentrations in the abutment screw and the abutment body.⁸⁵ Excessive occlusal stresses directed at the abutment fixation screw have been cited as a leading cause of screw loosening.⁸⁶ If left untreated, loose abutment screws can lead to crestal bone loss that would further expose the roughened surface to the biologic width.

In an animal model, another study evaluated the influence of the implant-abutment microgap on periimplant tissues in the mandibles of 4 canines.⁸⁷ The researchers placed implants and abutments in matching diameters 1 mm above, 1 mm below (countersunk), and level with the crestal bone.⁸⁷ After a 3-month nonsubmerged healing period, abutments were placed on the implants and the animals were allowed to function for 3 months before sacrifice and histologic evaluation.⁸⁷ Results showed that not only did placing the implant-abutment microgap deeper in the bone not result in additional bone loss, but that implants with countersunk microgaps had the least bone loss of the 3 study groups, which challenged the microgap theory.⁸⁷ Results of comparative short-term clinical studies have been inconclusive regarding the clinical efficacy of platform switching.^{88, 89} Additional long-term comparative data are still needed before definitive conclusions can be drawn to adequately support evidence-based treatment planning. Thus, while researchers have proposed platform shifting and surgical positioning of the implant collar at or above the bony crest, there is little data to suggest that these strategies will prevent all crestal bone loss around the implant. Therefore, maintaining 1.0 mm machined surface on the implant's cervical collar can help to resist the risk of bacterial attachment in the event of exposure to the gingival crevice or the oral cavity itself.

PERIIMPLANTITIS RISKS

Bacteria are significantly involved in the etiology of periimplant diseases, which can range from implant peri-mucositis to full periimplantitis that affects the supporting tissues around a dental implant and abutment.³⁰ The periimplantitis process typically begins as mucositis after 10-14 days of plaque retention, and triggers the same pathogen-induced inflammatory progression as periodontitis.³⁰ Consequently, patients with a history of periodontitis remain periodontitis-susceptible and have a higher risk of developing periimplantitis and implant failure than patients with no history of the disease.^{30, 90, 91} An estimated 16-28% of all implant patients will develop peri-implantitis lesions with marginal bone loss of at least 2 mm after 5-10 years.⁹⁰ This same disease process will also occur at a higher prevalence among patients with multiple implants than with an implant-supported, single tooth restoration.⁹⁰ For long-term implant health, and especially in periodontitis-susceptible patients, maintaining a 1.0 mm machined collar on dental implants may thus help to avoid roughened surface exposure and help to mitigate bacterial plaque adhesion to the dental implants, especially in periodontitis-susceptible patients.

REFERENCES

1. Pjetursson B, Asgeirsson AG, Zwahlen M, Sailer I. Improvements in implant dentistry over the last decade: comparison of survival and complication rates in older and newer publications. *Int J Oral Maxillofac Implants.* 2014;29(Suppl):308-324.
2. Glauser R. Implants with an oxidized surface placed predominantly in soft bone quality and subjected to immediate occlusal loading: results from an 11-year clinical follow-up. *Clin Implant Dent Relat Res.* 2015 Aug 13. doi: 10.1111/cid.12327. [Epub ahead of print].
3. Vandeweghe S, Hawker P, De Bruyn H. An up to 12-year retrospective follow-up on immediately loaded, surface-modified implants in the edentulous mandible. *Clin Implant Dent Relat Res.* 2015 Apr 23. doi: 10.1111/cid.12322. [Epub ahead of print].
4. Harel N, Piek D, Livne S, Palti A, Ormianer Z. A 10-year retrospective clinical evaluation of immediately loaded tapered maxillary implants. *Int J Prosthodont.* 2013;26:244-249.
5. Chen H, Liu N, Xu X, Qu X, Lu E. Smoking, radiotherapy, diabetes and osteoporosis as risk factors for dental implant failure: a meta-analysis. *PLoS ONE* 8(8): e71955. doi:10.1371/journal.pone.0071955. (URL at: <http://www.ncbi.nlm.nih.gov/pmc/articles/PMC3733795/>; last accessed May 4, 2016).
6. Hwang D, Wang H-L. Medical contraindications to implant therapy: part I: absolute contraindications. *Implant Dent.* 2006;15:353-360.
7. Hwang D, Wang H-L. Medical contraindications to implant therapy: part II: relative contraindications. *Implant Dent.* 2007;16:13-23.
8. Paquette DW, Brodala N, Williams RC. Risk factors for endosseous dental implant failure. *Dent Clin N Am.* 2006;50:361-374.
9. Sugerman PB, Barber MT. Patient selection for endosseous dental implants: oral and systemic considerations. *Int J Oral Maxillofac Implants.* 2002;17:191-201.
10. Chanavaz M. Patient screening and medical evaluation for implant and preprosthetic surgery. *J Oral Implantol.* 1998;24:222-229.
11. Lin S, Levin L, Goldman S, Peleg K. Dento-alveolar and maxillofacial injuries: a 5-year multi-center study. Part 1: general vs facial and dental trauma. *Dental Traumatology.* 2008;24:53-55.
12. Lin S, Levin L, Goldman S, Sela G. Dento-alveolar and maxillofacial injuries: a 5-year multi-center study. Part 2:

- severity and location. *Dental Traumatology*. 2008;24:56-58.
13. Giglio JA, Laskin DM. Perioperative errors contributing to implant failure. *Oral Maxillofac Surg Clin North Am*. 1998;10:197-202.
14. Olmedo D, Fernández MM, Guglielmotti MB, Cabrini RL. Macrophages related to dental implant failure. *Implant Dent*. 2003;12:75-80.
15. Camps-Font O, Figueiredo R, Valmaseda-Castellón E, Gay-Escoda C. Postoperative infections after dental implant placement: prevalence, clinical features, and treatment. *Implant Dent*. 2015;24:713-719.
16. Mazor Z, Cohen DK. Preliminary 3-dimensional surface texture measurement and early loading results with a microtextured implant surface. *Int J Oral Maxillofac Implants*. 2003;18:729-738.
17. Shalabi MM, Gortemaker A, Van't Hof MA, Jansen JA, Creugers NHJ. Implant surface roughness and bone healing: a systematic review. *J Dent Res*. 2006;85:496-500.
18. Trisi P, Lazzara R, Rao W, Rebaudi A. Bone-implant contact and bone quality: evaluation of expected and actual bone contact on machined and Osseotite implant surfaces. *Int J Periodontics Restorative Dent*. 2002;22:535-545.
19. Trisi P, Marcato C, Todisco M. Bone-to-implant apposition with machined and MTX microtextured implant surfaces in human sinus grafts. *Int J Periodontics Restorative Dent*. 2003;23:427-437.
20. Al-Nawas B, Götz H. Three-dimensional topographic and metrologic evaluation of dental implants by confocal laser scanning microscopy. *Clin Implant Dent Relat Res*. 2003;5:176-183.
21. Wennerberg A, Albrektsson T. On implant surfaces: a review of current knowledge and opinions. *Int J Oral Maxillofac Implants*. 2009;24:63-74.
22. Adell R. Long-term treatment results. Chapter 10. In: Brånemark P-I, Zarb GA, Albrektsson T (eds): *Tissue-Integrated Prostheses. Osseointegration in Clinical Dentistry*. Chicago: Quintessence Publishing Co., Inc., 1985:175-186.
23. Lindquist LW, Carlsson GE, Jemt T. A prospective 15-year follow-up study of mandibular fixed prostheses supported by osseointegrated implants. Clinical results and marginal bone loss. *Clin Oral Implants Res*. 1996;7:329-336.
24. Åstrand P, Engquist B, Dahlgren S, Gröndahl K, Engquist E, Feldmann H. Astra Tech and Brånemark system implants: a 5-year prospective study of marginal bone reactions. *Clin Oral Implants Res*. 2004;15:413-420.
25. Koirala DP, Singh SV, Chand P, Siddharth R, Jurel SK, Agarwal H, et al. Early loading of delayed versus immediately placed implants in the anterior mandible: a pilot comparative clinical study. *J Prosthet Dent*. 2016 Apr 22. pii: S0022-3913(16)00152-9. doi: 10.1016/j.prodent.2016.02.011. [Epub ahead of print].
26. Lang NP, Pun L, Lau KY, Li KY, Wong MCM. A systematic review on survival and success rates of implants placed immediately into fresh extraction sockets after at least 1 year. *Clin Oral Implants Res* 2015;23(Suppl 5):39-66.
27. Lee CT, Chiu TS, Chuang SK, Tarnow D, Stroupe J. Alterations of the bone dimension following immediate implant placement into extraction socket: systematic review and meta-analysis. *J Clin Periodontol* 2014;41:914-926.
28. Cohan FM. What are bacterial species? *Annu Rev Microbiol* 2002;56:457-458.
29. Siqueira JF Jr, Rocas IN. The oral microbiota: General overview, taxonomy, and nucleic acid techniques. *Methods Mol Biol* 2010;666:55-69.

30. Ormianer Z, Palti A. The use of tapered implant in the maxillae of periodontally susceptible patients: 10-year outcomes. *Int J Oral Maxillofac Implants* 2012;27:442-448.
31. Nunn ME. Understanding the etiology of periodontitis: An overview of periodontal risk factors. *Periodontology* 2000 2003;321:1-23.
32. Kaufman E. The new classification system of periodontal diseases and conditions. *Dent Today* 2001;20:102-105.
33. Lee A, Wang HL. Biofilm related to dental implants. *Implant Dent* 2010;19:387-393.
34. Oh TJ, Yoon J, Misch CE, Wang HL. The causes of early implant bone loss: myth or science? *J Periodontol*. 2002;73:322-333.
35. Dominici JT. Prosthodontic considerations in first stage implant failures. *Oral Maxillofac Surg Clin North Am* 1998;10:235-274.
36. Lewis SG. Prosthodontic considerations in implant failure. *Oral Maxillofac Surg Clin North Am*. 1998;10:309-322.
37. Ericksson RA, Albrektsson T. The effect of heat on bone regeneration. *J Oral Maxillofac Surg*. 1984;42:701-711.
38. Adell R, Lekholm U, Rockler B, Branemark P-I. A 15-year study of osseointegrated implants in the treatment of the edentulous jaw. *Int J Oral Surg*. 1981;10:387-416.
39. Misch CE, Dietsh-Misch F, Hoar J, Beck G, Hazen R, Misch CM. A bone quality-based implant system: first year of prosthetic loading. *J Oral Implantol*. 1999;25:185-197.
40. Satomi K, Akagawa Y, Nikai H, Tsuru H. Bone-implant interface structures after nontapping and tapping insertion of screw-type titanium alloy endosseous implants. *J Prosthet Dent*. 1988;59:339-342.
41. Haraldson T. A photoelastic study of some biomechanical factors affecting the anchorage of osseointegrated implants in the jaw. *Scand J Plast Reconstr Surg*. 1980;14:209-214.
42. Adell R, Lekholm U, Rockler B, et al. Marginal tissue reactions at osseointegrated titanium fixtures (I): a 3-year longitudinal prospective study. *Int J Oral Maxillofac Surg*. 1986;15:39-52.
43. Brocard D, Barthet P, Baysse E, et al. A multicenter report on 1,022 consecutively placed ITI implants: a 7-year longitudinal study. *Int J Oral Maxillofac Implants*. 2000;15:691-700.
44. Chrcanovic BR, Albrektsson T, Wennerberg A (2014). Flapless versus conventional flapped dental implant surgery: a meta-analysis. *PLoS ONE* 9(6): e100624. doi:10.1371/journal.pone.0100624. (URL at: <http://www.ncbi.nlm.nih.gov/pmc/articles/PMC4065043/pdf/pone.0100624.pdf>; last accessed on 5 May 2016).
45. Gotfredsen K. A 5-year prospective study of single-tooth replacements supported by the Astra Tech® implant: a pilot study. *Clin Implant Dent Relat Res*. 2004;6:1-8.
46. Schnitman PA, Shulman LB. Recommendations of the consensus development conference on dental implants. *J Am Dent Assoc*. 1979;98:373-377.
47. Albrektsson T, Isidor F. Consensus report of session V. In: Lang NP, Karring T, eds. *Proceedings of the 1st European Workshop on Periodontology*. London: Quintessence; 1993:365-369.
48. Zarb GA, Albrektsson T. Consensus report: towards optimized treatment outcomes for dental implants. *J Prosthet Dent*. 1998;80:641.

49. Quirynen M, Vogels R, Peeters W, van Steenberghe D, Naert I, Haffajee A. Dynamics of initial subgingival colonization of "pristine" peri-implant pockets. *Clin Oral Implants Res.* 2005;17:25-37.
50. Quirynen M, Teughels W. Microbiologically compromised patients and impact on oral implants. *Periodontology* 2000. 2003;33:119-128.
51. Sumida S, Ishihara K, Kishi M, Okuda K. Transmission of periodontal disease-associated bacteria from teeth to osseointegrated implant regions. *Int J Oral Maxillofac Implants.* 2002;17:696-702.
52. Heydenrijk K, Meijer HJA, van der Reijden WA, Raghoobar GM, Vissink A, Stengenga B. Microbiota around root-form endosseous implants: a review of the literature. *Int J Oral Maxillofac Implants.* 2002;17:829-838.
53. Covani U, Marconcini S, Crespi R, Barone A. Bacterial plaque colonization around dental implant surfaces. *Implant Dent.* 2006;15:298-304.
54. Cochran DL. Inflammation and bone loss in periodontal disease. *J Periodontol.* 2008;79:1569-1576.
55. Graves D. Cytokines that promote periodontal tissue destruction. *J Periodontol.* 2008;79:1585-1591.
56. Heckman SM, Linke JJ, Graef F, Foitzik Ch, Wichmann MG, Weber H-P. Stress and inflammation as a detrimental combination for peri-implant bone loss. *J Dent Res.* 2006;85:711-716.
57. Ferreira CF, Buttendorf AR, de Souza JG, Dalago H, Guenther SF, Bianchini MA. Prevalence of peri-implant diseases: analyses of associated factors. *Eur J Prosthodont Restor Dent.* 2015;23:199-206.
58. Misch CE, Suzuki JB, Misch-Dietsh FM, Bidez MW. A positive correlation between occlusal trauma and peri-implant bone loss: literature support. *Implant Dent.* 2005;14:108-116.
59. De Smet E, van Steenberghe D, Quirynen M, Naert I. The influence of plaque and/or excessive loading on marginal soft and hard tissue reactions around Branemark implants: a review of literature and experience. *Int J Periodontics Restorative Dent.* 2001;21:381-393.
60. Miyata T, Kobayashi Y, Araki H, Ohto T, Shin K. The influence of controlled occlusal overload on peri-implant tissue. Part 3: a histologic study in monkeys. *Int J Oral Maxillofac Implants.* 2000;15:425-431.
61. Saadoun AP, Le Gall M, Kricheck M. Microbial infections and occlusal overload: causes of failure in osseointegrated implants. *Pract Periodontics Aesthet Dent.* 1993;5:11-20.
62. Spray JR, Black CG, Morris HF, et al. The influence of bone thickness on facial marginal bone response: stage 1 placement through stage 2 uncovering. *Ann Periodontol.* 2000;5:119-128.
63. Penarrocha M, Palomar M, Sanchis SM, Guarinos J, Balaguer J. Radiographic study of marginal bone loss around 108 dental implants and its relationship to smoking, implant location, and morphology. *Int J Oral Maxillofac Implants.* 2004;19:861-867.
64. Galindo-Moreno P, Fauri M, Avila-Ortiz G, Fernandez- Barbero JE, Cabrera-Leon A, Sanchez-Fernandez E. Influence of alcohol and tobacco habits on peri-implant marginal bone loss: a prospective study. *Clin Oral Implants Res.* 2005;16:579-587.
65. Truhlar RS, Orenstein IH, Morris HF, Ochi S. Distribution of bone quality in patients receiving endosseous dental implants. *J Oral Maxillofac Surg.* 1997;55(suppl 5):38-45.
66. Petrie CS, Williams JL. Probabilistic analysis of peri-implant strain predictions as influenced by uncertainties in bone properties and occlusal forces. *Clin Oral Implants Res.* 2007;18:611-619.

67. Sennerby L, Gottlow J. Clinical outcomes of immediate/early loading of dental implants. A literature review of recent controlled prospective clinical studies. *Aust Dent J.* 2008;53(Suppl 1):S82-S88.
68. Fazel A, Aalai S, Rismanchian M. Effect of macro-design of immediately loaded implants on micromotion and stress distribution in surrounding bone using finite element analysis. *Implant Dent.* 2009;18:345-352.
69. Winter W, Klein D, Karl M. Effect of model parameters on finite element analysis of micromotions in implant dentistry. *J Oral Implantol.* 2013;39:23-29.
70. Brunette DM, Chehroudi B. The effects of the surface topography of micromachined titanium substrata on cell behavior in vitro and in vivo. *J Biomech Eng.* 1999 Feb;121(1):49-57.
71. Chehroudi B, Gould TRL, Brunette DM. Effects of a grooved titanium-coated implant surface on epithelial cell behavior in vitro and in vivo. *J Biomed Mater Res.* 1989;23:1067-1085.
72. Lee SW, Lee MH, Oh N, Park JA, Leesungbok R, Ahn SJ. Hydrophilicity and osteoblastic differentiation on microgrooved titanium substrata. *J Oral Implantol.* 2012;38:11-19.
73. Hermann JS, Jones AA, Bakaeen LG, Buser D, Schoolfield JD, Cochran DL. Influence of a machined collar on crestal bone changes around titanium implants: a histometric study in the canine mandible. *J Periodontol.* 2011;82:1329-1338.
74. Shin SY, Han DH. Influence of a microgrooved collar design on soft and hard tissue healing of immediate implantation in fresh extraction sites in dogs. *Clin Oral Implants Res.* 2010;21:804-814.
75. Bae HUK, Chung MK, Cha IH, Han DH. Marginal tissue response to different implant neck design. *J Korean Acad Prosthodont.* 2008;46:602-609.
76. Lee DW, Choi YS, Park KH, Kim CS, Moon IS. Effect of microthread on the maintenance of marginal bone level: a 3-year prospective study. *Clin Oral Implants Res.* 2007;18:465-470.
77. Nickenig HJ, Wichmann M, Schlegel KA, Nkenke E, Eitner S. Radiographic evaluation of marginal bone levels adjacent to parallel-screw cylinder machined-neck implants and rough-surfaced microthreaded implants using digitized panoramic radiographs. *Clin Oral Implants Res.* 2009;20:550-554.
78. Bratu EA, Tandlich M, Shapira L. A rough surface implant neck with microthreads reduces the amount of marginal bone loss: a prospective clinical study. *Clin Oral Implants Res.* 2009;20:827-832.
79. Chehroudi B, Gould TRL, Brunette DM. Titanium-coated micromachined grooves of different dimensions affect epithelial and connective-tissue cells differently in vivo. *J Biomed Mater Res.* 1990;24:1203-1219.
80. Abrahamsson I, Berglundh T. Tissue characteristics at microthreaded implants: an experimental study in dogs. *Clin Implant Dent Relat Res.* 2006;8:107-113.
81. Lee DW, Choi YS, Park KH, Kim CS, Moon IS. Effect of microthread on the maintenance of marginal bone level: a 3-year prospective study. *Clin Oral Implants Res.* 2007;18:465-470.
82. Astrand P, Engquist B, Dahlgren S, Grondahl K, Engquist E, Feldmann H. Astra Tech and Branemark system implants: a 5-year prospective study of marginal bone reactions. *Clin Oral Implants Res.* 2004;15:413-420.
83. Den Hartog, Meijer HJA, Tymstra N, Vissink A, Raghoobar GM. Single implants with different neck designs in the aesthetic zone: a randomized clinical trial. *Clin Oral Implants Res.* 2011;22:1289-1297.

84. Zitzmann NU, Abrahamsson I, Berglundh T, Lindhe J. Soft tissue reactions to plaque formation at implant abutments with different surface topography. An experimental study in dogs. *J Clin Periodontol.* 2002;29:456-461.
85. Maeda Y, Miura J, Taki I, Sogo M. Biomechanical analysis on platform switching: is there any biomechanical rationale? *Clin Oral Implants Res.* 2007;18:581-584.
86. Binon PP. Implants and components: entering the new millennium. *Int J Oral Maxillofac Implants.* 2000;15:76-94.
87. Todescan FF, Pustiglioni FE, Imbronito AV, Albrektsson T, Gioso M. Influence of the microgap in the peri-implant hard and soft tissues: a histomorphometric study in dogs. *Int J Oral Maxillofac Implants.* 2002;15:467-472.
88. Vela-Nebot X, Rodriguez-Ciurana X, Rodado-Alonso C, Segala-Torres M. Benefits of an implant platform modification technique to reduce crestal bone resorption. *Implant Dent.* 2006;15:313-320.
89. Becker J, Ferrari D, Herten M, Kirsch A, Schaer A, Schwarz F. Influence of platform switching on crestal bone changes at non- submerged titanium implants: a histomorphometrical study in dogs. *J Clin Periodontol.* 2007;34:1089-1096.
90. Lundgren D, Rylander H, Laurell L. To save or to extract, that is the question. Natural teeth or dental implants in periodontitis-susceptible patients: clinical decision-making and treatment strategies exemplified with patient case presentations. *Periodontology 2000.* 2008;47:27-50.
91. Van der Weijden GA, van Bommel KM, Renvert S. Implant therapy in partially edentulous periodontally compromised patients: a review. *J Clin Periodontol.* 2005;32:506-511.

TABLES

Source	Type	Examples
Manufacturer	Manufacturing errors	Metal burs from machining Surface contamination (e.g. manufacturing oil, particulate debris) Wrong labeling
	Design errors	Interfacial gaps between machined components Abutment rotational instability Inadequate or poor engagement of screw threads
Patient	Uncontrolled diseases	Periodontitis Diabetes mellitus Osteoporosis
	Compromising medications	Intravenous bisphosphonates Long-term corticosteroid use
	Parafunctional habits	Bruxing Fingernail biting / pencil chewing
	Compromising personal habits	Smoking Alcoholism Poor oral hygiene
Clinician	Latrogenic errors	Excessive drilling heat Improper implant angulation Ineffective occlusal scheme

*Examples are not an exhaustive list; other examples also apply

Time Interval	1979 [†]	1985 [‡]	1986 [‡]	1993 [‡]	1998 [‡]
General	1/3 of implant height ⁴⁶	Not Reported	Not Reported	Not Reported	Not Reported
After submerged healing	Not Reported	Not Reported	Not Reported	Not Reported	Not Reported
First year of loading	Not Reported	1.0-1.5 mm ²²	<1.0 mm ⁴²	<1.0 mm ⁴⁷	Not Reported
Follow-up years	Not Reported	Not Reported	0.05-0.01 ⁴²	<0.2 mm ⁴⁷	<0.2 mm ⁴⁸

[†]Blade implants only

[‡]Root-form implants only